

# DDX51 Antibody



PACO08796

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DDX51

**Uniprot**

Q8N8A6

**Synonyms:**

DEAD (Asp-Glu-Ala-Asp) box polypeptide 51; DDX51; DKFZp686N2081; MGC42193 ;

**Immunogen:**

Human DDX51.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A