

DEDD Antibody



PACO08800

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DEDD

Uniprot

O75618

Synonyms:

death effector domain containing; DEDD; CASP8IP1; DEDD1; DEFT; FLDED1; KE05 ;

Immunogen:

Human DEDD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A