

# DFNA5 Antibody



PACO08808

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GSDME

**Uniprot**

O60443

**Synonyms:**

deafness, autosomal dominant 5; DFNA5; ICERE-1 ;

**Immunogen:**

Human DFNA5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A