

DHRS11 Antibody



PACO08824

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DHRS11

Uniprot

Q6UWP2

Synonyms:

dehydrogenase/reductase (SDR family) member 11; DHRS11; ARPG836; FLJ39232; MGC4172; SDR24C1 ;

Immunogen:

Human DHRS11.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A