

DHX15 Antibody



PACO08829

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DHX15

Uniprot

O43143

Synonyms:

DEAH (Asp-Glu-Ala-His) box polypeptide 15; DHX15; DBP1; DDX15; HRH2; PRP43; PRPF43; PrPp43p ;

Immunogen:

Human DHX15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A