

DHX16 Antibody



PACO08830

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DHX16

Uniprot

O60231

Synonyms:

DEAH (Asp-Glu-Ala-His) box polypeptide 16; DHX16; DBP2; DDX16; PRO2014; PRP8 ;

Immunogen:

Human DHX16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A