

DHX38 Antibody



PACO08834

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DHX38

Uniprot

Q92620

Synonyms:

DEAH (Asp-Glu-Ala-His) box polypeptide 38; DHX38; DDX38; KIAA0224; PRP16; PRPF16 ;

Immunogen:

Human DHX38.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A