

DIDO1 Antibody



PACO08840

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DIDO1

Uniprot

Q9BTC0

Synonyms:

death inducer-obliterator 1; DIDO1; BYE1; C20orf158; DATF1; DIDO2; DIDO3; DIO-1; DIO1; DKFZp434P1115; FLJ11265; KIAA0333; MGC16140; dJ885L7.8 ;

Immunogen:

Human DIDO1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A