

DLAT Antibody



PACO08848

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DLAT

Uniprot

P10515

Synonyms:

dihydrolipoamide S-acetyltransferase; DLAT; DLTA; PDC-E2; PDCE2 ;

Immunogen:

Human DLAT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A