

DLG3 Antibody



PACO08853

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DLG3

Uniprot

Q92796

Synonyms:

discs, large homolog 3 (Drosophila); DLG3; KIAA1232; MRX; MRX90; NE-Dlg; NEDLG; SAP-102; SAP102 ;

Immunogen:

Human DLG3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A