

# DMC1 Antibody



PACO08863

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DMC1

**Uniprot**

Q14565

**Synonyms:**

DMC1 dosage suppressor of mck1 homolog, meiosis-specific homologous recombination (yeast); DMC1; DMC1H; HsLim15; LIM15; MGC150472; MGC150473; dJ199H16.1 ;

**Immunogen:**

Human DMC1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A