

DMWD Antibody



PACO08866

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DMWD

Uniprot

Q09019

Synonyms:

dystrophia myotonica, WD repeat containing; DMWD; D19S593E; DMR-N9; DMRN9; gene59 ;

Immunogen:

Human DMWD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A