

DNAJB1 Antibody



PACO08874

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DNAJB1

Uniprot

P25685

Synonyms:

DnaJ (Hsp40) homolog, subfamily B, member 1; DNAJB1; HSPF1; Hdj1; Hsp40; Sis1 ;

Immunogen:

Human DNAJB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A