

DNAJB14 Antibody



PACO08877

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNAJB14

Uniprot

Q8TBM8

Synonyms:

DnaJ (Hsp40) homolog, subfamily B, member 14; DNAJB14; EGNR9427; FLJ14281; MGC22187; PRO34683 ;

Immunogen:

Human DNAJB14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A