

DNAJB3 Antibody



PACO08879

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DNAJB3

Uniprot

Q8WWF6

Synonyms:

DnaJ (Hsp40) homolog, subfamily B, member 3; DNAJB3; HCG3; MGC26879 ;

Immunogen:

Human DNAJB3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A