

DNAJC15 Antibody



PACO08887

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DNAJC15

Uniprot

Q9Y5T4

Synonyms:

DnaJ (Hsp40) homolog, subfamily C, member 15; DNAJC15; DNAJD1; HSD18; MCJ ;

Immunogen:

Human DNAJC15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A