

# DNAJC2 Antibody



PACO08890

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DNAJC2

**Uniprot**

Q99543

**Synonyms:**

DnaJ (Hsp40) homolog, subfamily C, member 2; DNAJC2; MPHOSPH11; MPP11; ZRF1; ZUO1 ;

**Immunogen:**

Human DNAJC2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A