

DNASE2 Antibody



PACO08900

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DNASE2

Uniprot

O00115

Synonyms:

deoxyribonuclease II, lysosomal; DNASE2; DNASE2A; DNL; DNL2 ;

Immunogen:

Human DNASE2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A