

DOK6 Antibody



PACO08917

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DOK6

Uniprot

Q6PKX4

Synonyms:

docking protein 6; DOK6; DOK5L; HsT3226; MGC20785 ;

Immunogen:

Human DOK6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A