

PACO08918

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DXO

Uniprot

O77932

Synonyms:

dom-3 homolog Z (C. elegans); DOM3Z; DOM3L; NG6 ;

Immunogen:

Human DOM3Z.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A