

DPF2 Antibody



PACO08920

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

DPF2

Uniprot

Q92785

Synonyms:

D4, zinc and double PHD fingers family 2; DPF2; MGC10180; REQ; UBID4; ubi-d4 ;

Immunogen:

Human DPF2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A