

# DTD1 Antibody



PACO08942

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DTD1

**Uniprot**

Q8TEA8

**Synonyms:**

D-tyrosyl-tRNA deacylase 1 homolog (S. cerevisiae); DTD1; C20orf88; DUEB; HARS2; MGC119131; MGC41905; bA379J5.3; bA555E18.1; pqn-68 ;

**Immunogen:**

Human DTD1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A