

DTNA Antibody



PACO08944

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DTNA

Uniprot

Q9Y4J8

Synonyms:

dystrobrevin, alpha; DTNA; D18S892E; DRP3; DTN; FLJ96209; LVNC1 ;

Immunogen:

Human DTNA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A