

DTX2 Antibody



PACO08946

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DTX2

Uniprot

Q86UW9

Synonyms:

deltex homolog 2 (Drosophila); DTX2; KIAA1528; MGC71098; RNF58 ;

Immunogen:

Human DTX2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A