

DUSP11 Antibody



PACO08952

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DUSP11

Uniprot

O75319

Synonyms:

dual specificity phosphatase 11 (RNA/RNP complex 1-interacting); DUSP11; PIR1 ;

Immunogen:

Human DUSP11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A