

# DUSP13 Antibody



PACO08953

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DUSP13

**Uniprot**

Q6B8I1

**Synonyms:**

dual specificity phosphatase 13; DUSP13; BEDP; DUSP13A; DUSP13B; FLJ32450; MDSP; SKRP4; TMDP ;

**Immunogen:**

Human DUSP13.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A