

DYNC2LI1 Antibody



PACO08965

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DYNC2LI1

Uniprot

Q8TCX1

Synonyms:

dynein, cytoplasmic 2, light intermediate chain 1; DYNC2LI1; CGI-60; D2LIC; DKFZp564A033; LIC3 ;

Immunogen:

Human DYNC2LI1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A