

E2F8 Antibody



PACO08976

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

E2F8

Uniprot

A0AVK6

Synonyms:

E2F transcription factor 8; E2F8; FLJ23311 ;

Immunogen:

Human E2F8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A