

EBF3 Antibody



PACO08980

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EBF3

Uniprot

Q9H4W6

Synonyms:

early B-cell factor 3; EBF3; COE3; O/E-2 ;

Immunogen:

Human EBF3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A