

EBP Antibody



PACO08982

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

EBP

Uniprot

Q15125

Synonyms:

emopamil binding protein (sterol isomerase); EBP; CDPX2; CHO2; CPX; CPXD ;

Immunogen:

Human EBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A