

ECD Antibody



PACO08983

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ECD

Uniprot

O95905

Synonyms:

ecdysoneless homolog (Drosophila); ECD; GCR2; HSGT1 ;

Immunogen:

Human ECD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A