

ECH1 Antibody



PACO08985

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ECH1

Uniprot

Q13011

Synonyms:

enoyl Coenzyme A hydratase 1, peroxisomal; ECH1; HPXEL ;

Immunogen:

Human ECH1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A