

EDAR Antibody



PACO08989

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EDAR

Uniprot

Q9UNE0

Synonyms:

ectodysplasin A receptor; EDAR; DL; ED1R; ED3; ED5; EDA-A1R; EDA1R; EDA3; FLJ94390; HRM1 ;

Immunogen:

Human EDAR.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A