

EEF1E1 Antibody



PACO09000

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EEF1E1

Uniprot

O43324

Synonyms:

eukaryotic translation elongation factor 1 epsilon 1; EEF1E1; AIMP3; P18 ;

Immunogen:

Human EEF1E1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A