

EGFLAM Antibody



PACO09014

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EGFLAM

Uniprot

Q63HQ2

Synonyms:

EGF-like, fibronectin type III and laminin G domains; EGFLAM; AGRINL; AGRNL; FLJ39155 ;

Immunogen:

Human EGFLAM.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A