

# EIF2AK2 Antibody



PACO09032

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EIF2AK2

**Uniprot**

P19525

**Synonyms:**

eukaryotic translation initiation factor 2-alpha kinase 2; EIF2AK2; EIF2AK1; MGC126524; PKR; PRKR ;

**Immunogen:**

Human EIF2AK2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A