

EIF2AK3 Antibody



PACO09033

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

EIF2AK3

Uniprot

Q9NZJ5

Synonyms:

eukaryotic translation initiation factor 2-alpha kinase 3; EIF2AK3; DKFZp781H1925; HRI; PEK; PERK; WRS ;

Immunogen:

Human PERK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A