

# AGO2 Antibody



PACO09039

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

AGO2

**Uniprot**

Q9UKV8

**Synonyms:**

eukaryotic translation initiation factor 2C, 2; EIF2C2; AGO2; MGC3183; Q10 ;

**Immunogen:**

Human EIF2C2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A