

EIF3M Antibody



PACO09050

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

EIF3M

Uniprot

Q7L2H7

Synonyms:

eukaryotic translation initiation factor 3, subunit M; EIF3M; B5; FLJ29030; GA17; PCID1; hfl-B5 ;

Immunogen:

Human EIF3M.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A