

# EIF4E3 Antibody



PACO09055

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

EIF4E3

**Uniprot**

Q8N5X7

**Synonyms:**

eukaryotic translation initiation factor 4E family member 3; EIF4E3; MGC39820; MGC86971 ;

**Immunogen:**

Human EIF4E3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A