

EIF4H Antibody



PACO09060

Product Information

Size:

50ul

Protein Background:**Reactivity:**

Human, Mouse, Rat

Gene ID:

EIF4H

Source:

Rabbit

Uniprot

Q15056

Isotype:

IgG

Synonyms:

27kd;

Applications:

ELISA

Immunogen:

Human EIF4H.

Recommended dilutions:**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A