

ENO1 Antibody



PACO09090

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ENO1

Uniprot

P06733

Synonyms:

enolase 1, (alpha); ENO1; ENO1L1; MBP-1; MPB1; NNE; PPH ;

Immunogen:

Human ENO1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A