

ENPP1 Antibody



PACO09096

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ENPP1

Uniprot

P22413

Synonyms:

ectonucleotide pyrophosphatase/phosphodiesterase 1; ENPP1; M6S1; NPP1; NPPS; PC-1; PCA1; PDNP1 ;

Immunogen:

Human ENPP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A