

PACO09100

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ENSA

Uniprot

O43768

Synonyms:

endosulfine alpha; ENSA; MGC4319; MGC78563; MGC8394 ;

Immunogen:

Human ENSA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A