

ENTPD1 Antibody



PACO09101

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ENTPD1

Uniprot

P49961

Synonyms:

ectonucleoside triphosphate diphosphohydrolase 1; ENTPD1; ATPDase; CD39; DKFZp686D194; DKFZp686I093; FLJ40921; FLJ40959; NTPDase-1 ;

Immunogen:

Human ENTPD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A