

EPB41 Antibody



PACO09106

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

EPB41

Uniprot

P11171

Synonyms:

erythrocyte membrane protein band 4.1 (elliptocytosis 1, RH-linked); EPB41;4.1R; EL1; HE ;

Immunogen:

Human EPB41.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A