

EPHX2 Antibody



PACO09113

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

EPHX2

Uniprot

P34913

Synonyms:

epoxide hydrolase 2, cytoplasmic; EPHX2; CEH; SEH ;

Immunogen:

Human EPHX2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A