

EPS15L1 Antibody



PACO09115

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

EPS15L1

Uniprot

Q9UBC2

Synonyms:

epidermal growth factor receptor pathway substrate 15-like 1; EPS15L1; EPS15R ;

Immunogen:

Human EPS15L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A