

EPS8L2 Antibody



PACO09117

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

EPS8L2

Uniprot

Q9H6S3

Synonyms:

EPS8-like 2; EPS8L2; EPS8R2; FLJ16738; FLJ21935; FLJ22171; MGC126530; MGC3088 ;

Immunogen:

Human EPS8L2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A