

ERAP1 Antibody



PACO09120

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ERAP1

Uniprot

Q9NZ08

Synonyms:

endoplasmic reticulum aminopeptidase 1; ERAP1; A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; KIAA0525; PILS-AP; PILSAP ;

Immunogen:

Human ERAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A